


Kaveh Moradkhani

Ph.D. student in Software and IT
Engineering at ÉTS University



In the world of medical image processing, our goal is to provide life-saving solutions and help those in need. It's a noble cause that we're proud to be a part of

Paper Presentation

October 18th 2023



Le génie pour l'industrie


Under the supervision of

Pro. Sylvain Bouix



Neuro-iX


Vox2Cortex: Fast Explicit Reconstruction of Cortical Surfaces from 3D MRI Scans with Geometric Deep Neural Networks



MAIN AIM



The main objective of this work is to develop a fast and highly accurate reconstruction of the cortex, which has high clinical value





UNDERSTANDING TOPIC



Unlocking the Power of Cortical Surface Reconstruction




- ❖ Precise analysis of cortical thickness
- ❖ Sulcal morphology




MAIN CHALLENGE








MAIN CHALLENGE



- ❖ Lengthy runtimes of multiple hours
- ❖ Intricate post-processing



FreeSurfer



RELATED WORK



- **Voxel-based surface reconstruction:**




- ❖ FastSurfer
- ❖ Deep implicit representations
- ❖ SegRecon

- **Mesh-based surface reconstruction:**

- ❖ Voxel2Mesh
- ❖ MeshDeformNet
- ❖ PIALNN



Related Work




METHODOLOGY



Detail of the individual building blocks:

- ❖ Image encoding and decoding (CNN)
- ❖ Mesh deformation (GNN)



❖ The loss function of Vox2Cortex:

$$\mathcal{L}(y^p, y^{gt}) = \mathcal{L}_{\text{vox}}(y^p, y^{gt}) + \mathcal{L}_{\text{mesh}}(y^p, y^{gt})$$

❖ Voxel loss:

$$\mathcal{L}_{\text{vox}}(y^p, y^{gt}) = \sum_{l=1}^L \mathcal{L}_{\text{BCE}}(B_l^p, B_l^{gt})$$



Mesh loss:

$$\mathcal{L}_{\text{mesh}}(y^{\text{P}}, y^{\text{gt}}) = \mathcal{L}_{\text{mesh}, \text{cons}}(y^{\text{P}}, y^{\text{gt}}) + \mathcal{L}_{\text{mesh}, \text{reg}}(y^{\text{P}}).$$

$$\begin{aligned} \mathcal{L}_{\text{mesh}, \text{cons}}(y^{\text{P}}, y^{\text{gt}}) = & \sum_{s=1}^S \sum_{c=1}^C \left[\lambda_{1,c} \mathcal{L}_{\text{C}}(\mathcal{M}_{s,c}^{\text{P}}, \mathcal{M}_c^{\text{gt}}) \right. \\ & \left. + \lambda_{2,c} \mathcal{L}_{\text{n}, \text{inter}}(\mathcal{M}_{s,c}^{\text{P}}, \mathcal{M}_c^{\text{gt}}) \right]. \end{aligned}$$

$$\begin{aligned} \mathcal{L}_{\text{mesh}, \text{reg}}(y^{\text{P}}) = & \sum_{s=1}^S \sum_{c=1}^C \left[\lambda_{3,c} \mathcal{L}_{\text{Lap}, \text{rel}}(\mathcal{M}_{s,c}^{\text{P}}, \Delta_{s,c}^{\text{P}}) \right. \\ & + \lambda_{4,c} \mathcal{L}_{\text{n}, \text{intra}}(\mathcal{M}_{s,c}^{\text{P}}) \\ & \left. + \lambda_{5,c} \mathcal{L}_{\text{edge}}(\mathcal{M}_{s,c}^{\text{P}}) \right]. \end{aligned}$$



DATASETS AND PRE-PROCESSING



❖ Pre-processing

- Registering MRI scans to the MNI152 space
- Resize to $128 \times 144 \times 128$ voxels


❖ Dataset

ADNI

- Provides MRI T1 scans for subjects with Alzheimer's Disease, Mild Cognitive Impairment, and healthy subjects
- Balanced according to diagnosis, age, and sex
- ADNIlarge contains 1,155 subjects
- ADNIsmall contains 299 subjects

OASIS

- Contains MRI T1 scans of 416 subjects
- 100 subjects have been diagnosed with very mild to moderate Alzheimer's disease
- Balanced on diagnosis, age, and sex



COMPARISON



- ❖ Average Symmetric Surface Distance (ASSD)
- ❖ 90-Percentile Hausdorff Distance (HD)

Data	Method	Left WM Surface		Right WM Surface		Left Pial Surface		Right Pial Surface	
		ASSD	HD	ASSD	HD	ASSD	HD	ASSD	HD
ADNI large	Vox2Cortex	0.345 ±0.056	0.720 ±0.125	0.347 ±0.046	0.720 ±0.087	0.327 ±0.031	0.755 ±0.102	0.318 ±0.029	0.781 ±0.102
	DeepCSR [4]	0.422 ±0.058	0.852 ±0.134	0.420 ±0.058	0.880 ±0.156	0.454 ±0.059	0.927 ±0.243	0.422 ±0.053	0.890 ±0.197
	nnUNet [23]	1.176 ±0.345	1.801 ±2.835	1.159 ±0.242	1.739 ±1.880	1.310 ±0.292	3.152 ±2.374	1.317 ±0.312	3.295 ±2.387
OASIS	Vox2Cortex	0.315 ±0.039	0.680 ±0.137	0.318 ±0.048	0.682 ±0.151	0.362 ±0.036	0.894 ±0.141	0.373 ±0.041	0.916 ±0.137
	DeepCSR [4]	0.360 ±0.042	0.731 ±0.104	0.335 ±0.050	0.670 ±0.195	0.458 ±0.056	1.044 ±0.290	0.442 ±0.058	1.037 ±0.294



LIMITATION



- ❖ No Ground Truth: they used FreeSurfer surfaces as pseudo ground-truth due to the absence of manually generated ground-truth surfaces for their data
- ❖ Limited Data Scope: their analysis only includes healthy subjects and those with dementia, so the model may not perform accurately in the presence of other brain morphological changes like tumors
- ❖ Data Specificity: their model's performance is based on the discussed data and may not generalize well to unseen data from different domains
- ❖ Fairness: they balanced data splits with respect to age and sex, their model might still have biases and potentially discriminate against underrepresented groups




FEATURE



- ❖ Vox2Cortex simultaneously predicts white matter (WM) and pial surfaces of both hemispheres, resulting in an output of four meshes
- ❖ Vox2Cortex models the interdependency between WM and pial surfaces by exchanging information between them
- ❖ In contrast to time-consuming traditional methods for brain analysis, such a network performs the segmentation in seconds
- ❖ The current standard for cortical surface reconstruction is FreeSufer, which produces smooth, accurate, and topologically correct surface meshes but runs for several hours per scan
- ❖ Current deep-learning approaches for cortical surface reconstruction focus on voxel-based or implicit surface reconstruction methods



Example



CONCLUSION



In summary, Vox2Cortex is a groundbreaking approach for explicit cortical surface reconstruction, offering speed, accuracy, and many promising applications



Thank you!

