

3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

Adnan Ahmad, Drew Parker, Suhani Dheer, Zahra Riahi Samani, Ragini Verma

**Computerized Medical Imaging and Graphics
2023**

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

- Quality control (QC) is crucial for artifact detection in diffusion MRI before any analysis
- Manual QC is time-consuming and subjective
- Need for automated QC pipeline
- Goal of the paper:


Proposing an automated deep learning pipeline for artifact detection in dMRI volumes



Using 3D-DenseNet across diverse dMRI datasets

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images


▪ Model architecture



- 3D-DenseNet with 3 dense blocks
- Transition layers for downsampling
- Global Average Pooling + softmax classifier
- Input volumes resized to 96×96×70

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images


■ Training details:



- Batch size: 5
- epochs: 20
- Optimizer: Adam
- Cross entropy loss
- Training time ~3h on NVIDIA 1080 GPU

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

▪ Inference pipeline:




- Accepting dMRI scans, preprocessing each 3D volume to a fixed size
- Running in inference mode, providing the user with predicted artifact / normal label for each volume, calculating and reporting performance metrics
- Generating a QC report listing flagged volumes so users can drop them before further analyses

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

▪ Dataset

➤ 7 datasets: 3 for training / validation, 4 for testing

➤ Total 9258 volumes from 678 subjects 

5619 training volumes
1292 validation volumes
2347 test volumes

➤ Heterogeneous in scanners and protocols

➤ Annotating by an expert with 8 years of experience identifying artifacts in diffusion scans

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

▪ Results

Results – 3D-QCNet Model.

	Dataset	Accuracy	Precision	Recall
Validation Set	Dataset 1	97	95	80
	Dataset 2	81	91	66
	Dataset 3	89	90	82
	<i>Average</i>	89	92	76
	Dataset 4	97	84	81
	Dataset 5	92	86	98
	Dataset 6	96	81	100
Test Set	Dataset 7	92	95	89
	<i>Average</i>	94	87	92

Defining “artifacts” as the positive class

✓ Consistent across 4 unseen datasets

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

■ Results




Fig. 3. Scans from the test set illustrated to demonstrate 3D-QCNet's model performance with respect to ground-truth. True Positive samples – A Ghosting artifact, B Herringbone artifact, C Motion/interslice instability artifact, D Faint Chemical artifact (marked in yellow). True Negative samples – E Weighted Image is noisy but is correctly marked as normal. F B0 image with no artifacts. False Positive samples – G Abnormal anatomy of the brain may be affecting the classifier. H Weighted image is noisy but there are no visible artifacts; the model may be too sensitive. False Negative samples – I Chemical shift artifact alongside instability and susceptibility.

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

▪ Results



Fig. 4. These scans are from a volume that was marked as normal by our annotator but 3D-QCNet labelled it as having an artifact. Later, on closer inspection it was found to have ghosting artifacts in the ventricles and subcortical regions along with some interslice instability.

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

- Results

Results on Dataset 5 when using labels from different annotators.

Ground Truth	Comparison	Accuracy	Precision	Recall
8 + years of experience ← Annotator 1	3D-QCNet	92.75	86.5	98.85
A month’s training on identifying artifacts ← Annotator 2	3D-QCNet	76	75.6	75.6

- ✓ The inherent inaccuracy of Annotator 2 given their relative inexperience
- ✓ Due to the data it was trained on, the model is more akin to Annotator 1 and learns their expertise as well as biases

❖ 3D-QCNet – A pipeline for automated artifact detection in diffusion MRI images

▪ Conclusion

- ✓ End-to-end automation of dMRI QC
- ✓ Robust generalization across diverse data
- × Binary output only. It flags “artifact” vs “normal” but doesn’t identify artifact type
- × Fixed input size (96×96×70)
- × Not training on data annotated by multiple experts with similar experience

Thank you