

Towards user-centered interactive medical image segmentation in VR with an assistive AI agent

Pascal Spiegler¹, Arash Harirpoush¹ and Yiming Xiao¹

Department of Computer Science and Software Engineering, Concordia University, Montreal, Quebec, Canada.

Challenges in segmentation


Problem:

- Manual segmentation of 3D medical scans is time-consuming, labor-intensive, and prone to errors
- Fully automatic algorithms lack reliability without user feedback
- Existing tools do not seamlessly integrate AI and VR for interactive, patient-specific segmentation

Goal:

- Develop an immersive VR system with AI assistance that facilitates **efficient, accurate, and user-friendly** segmentation
- Investigate **optimal interaction modes** (controller, head, eye) for mask refinement in VR

Solution proposed : SAMIRA



Development of SAMIRA, a conversational AI agent in VR for medical image segmentation.

Related works

1) VR & Medical Training:

Early VR systems support clinical education but lack AI integration for segmentation (*Chheang et al., 2024*)

2) AI in Medical Imaging:

Foundation models like **BiomedParse** and **SAM2** aid in segmentation but are limited to 2D or desktop applications (*Zhao et al., 2025, Kirillov et al., 2024*)

3) Human-Computer Interaction in VR:

- Studies compare input modes: controller, head pointing, eye tracking (*Sidenmark et al., 2023, Luro et al., 2019, Schultheis et al., 2024*)
- Prior work shows controllers excel in accuracy; head and eye tracking are more natural but less precise

Gap:

No prior VR system integrates **AI-powered, interactive 3D segmentation** with **human-in-the-loop workflows**


SAMIRA

Fusion of AI & VR: Combines foundational radiological models with an immersive VR interface for real-time, human-in-the-loop segmentation

Key Components:

- **AI Models:** BiomedParse and SAM2 perform initial segmentation, refinement, and mask propagation.
- **VR Interface:** Visualizes slices, reference images, and 3D reconstructions; supports speech and input modes for users.

Workflow : Interaction modes



The AI-Assisted Segmentation System

1) BiomedParse:

- A **Transformer-based vision-language model** trained on 1.2 million medical images with reports.
- It can **detect, classify, and segment** 82 clinical concepts (like tumors) across different imaging modalities using natural language prompts.



2) SAM2 (Segment Anything Model 2):

- A **model** designed for **interactive segmentation** using user prompts like points or boxes.
- It produces segmentation masks, which outline objects in images, and can **propagate these masks** across similar images or video frames

3) RAG (Retrieval-Augmented Generation):


An AI framework that **searches a database** for relevant example images or slices, then **uses those** to generate more accurate, case-specific responses.

Workflow steps




RAG Pipeline

RAG Request 1: General Contextualization of Target



RAG Request 2: Patient-Specific Guidance




Results

Table 2 Comparison of interaction paradigms for segmentation refinement. Values are mean \pm standard deviation. The best score is in bold. NASA-TLX is out of 100.

Interaction Mode	3D Dice (%)	Time (s)	NASA-TLX	Composite Score
Controller	99.25 \pm 0.25	220.3 \pm 79.3	18.8 \pm 14.5	0.51 \pm 1.91
Head Pointing	99.21 \pm 0.30	248.8 \pm 78.5	16.8 \pm 13.9	0.20 \pm 1.56
Eye Tracking	99.13 \pm 0.46	251.1 \pm 78.7	26.6 \pm 15.4	-0.71 \pm 1.77

Results



Conclusion

Key Findings:

- User-centered VR interaction paradigms are feasible and effective for medical segmentation
- AI assistance enhances confidence, understanding, and accuracy.
- Head Pointing offers a low mental effort alternative, ideal for single-hand use.

Implications:

- Potential to improve clinical workflows, training, and patient-specific analysis.
- Combines immersive visualization with intelligent guidance for better decision-making.

Future works:

- Further refine interaction techniques for real-world clinical settings.
- Expand to other medical applications and datasets.
- Explore integration with emerging VR and AI technologies.

Thank you!